

AP1M2 Antibody



PACO07765

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AP1M2

Uniprot

Q9Y6Q5

Synonyms:

adaptor-related protein complex 1, mu 2 subunit; AP1M2; AP1-mu2; HSMU1B; MU-1B; MU1B; mu2 ;

Immunogen:

Human AP1M2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A