

AP2B1 Antibody



PACO07767

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

AP2B1

Uniprot

P63010

Synonyms:

adaptor-related protein complex 2, beta 1 subunit; AP2B1; ADTB2; AP105B; AP2-beta; CLAPB1; DKFZp781K0743 ;

Immunogen:

Human AP2B1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A