

AP3D1 Antibody



PACO07769

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AP3D1

Uniprot

O14617

Synonyms:

adaptor-related protein complex 3, delta 1 subunit; AP3D1; ADTD; hBLVR ;

Immunogen:

Human AP3D1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A