

# AP3M2 Antibody



PACO07771

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

AP3M2

**Uniprot**

P53677

**Synonyms:**

adaptor-related protein complex 3, mu 2 subunit; AP3M2; AP47B; CLA20; P47B ;

**Immunogen:**

Human AP3M2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A