

APBA2 Antibody



PACO07775

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

APBA2

Uniprot

Q99767

Synonyms:

amyloid beta (A4) precursor protein-binding, family A, member 2; APBA2; D15S1518E; HsT16821; LIN-10; MGC99508; MGC:14091; MINT2; X11L ;

Immunogen:

Human APBA2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A