

APLF Antibody



PACO07783

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

APLF

Uniprot

Q8IW19

Synonyms:

aprataxin and PNKP like factor; APLF; APFL; C2orf13; FLJ16593; MGC47799; PALF; Xip1 ;

Immunogen:

Human APLF.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A