

APOBEC4 Antibody



PACO07792

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

APOBEC4

Uniprot

Q8WW27

Synonyms:

apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative);
APOBEC4; C1orf169; MGC26594 ;

Immunogen:

Human APOBEC4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A