

APP Antibody



PACO07797

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

APP

Uniprot

P05067

Synonyms:

amyloid beta (A4) precursor protein; APP; AAA; Abeta; ABPP; AD1; APPI; CTFgamma; CVAP; PN2 ; ; amyloid P; Amyloid beta A4 protein; A4; beta Amyloid

Immunogen:

Human APP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A