

NAA10 Antibody



PACO07803

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NAA10

Uniprot

P41227

Synonyms:

ARD1 homolog A, N-acetyltransferase (S. cerevisiae); ARD1A; ARD1; DXS707; MGC71248; TE2 ;

Immunogen:

Human ARD1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A