

# ARFIP2 Antibody



PACO07813

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARFIP2

**Uniprot**

P53365

**Synonyms:**

ADP-ribosylation factor interacting protein 2; ARFIP2; POR1 ; Arfaptin-2

**Immunogen:**

Human ARFIP2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A