

# ARHGAP15 Antibody



PACO07819

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARHGAP15

**Uniprot**

Q53QZ3

**Synonyms:**

Rho GTPase activating protein 15; ARHGAP15; BM046 ;

**Immunogen:**

Human ARHGAP15.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A