

ARHGAP18 Antibody



PACO07820

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ARHGAP18

Uniprot

Q8N392

Synonyms:

Rho GTPase activating protein 18; ARHGAP18; FLJ25728; MGC126757; MGC138145; MacGAP; bA307O14.2 ;

Immunogen:

Human ARHGAP18.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A