

ARHGAP20 Antibody



PACO07821

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGAP20

Uniprot

Q9P2F6

Synonyms:

Rho GTPase activating protein 20; ARHGAP20; KIAA1391; RARHOGAP ;

Immunogen:

Human ARHGAP20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A