

ARHGAP24 Antibody



PACO07822

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Dog, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ARHGAP24

Uniprot

Q8N264

Synonyms:

Rho GTPase activating protein 24; ARHGAP24; DKFZp564B1162; FILGAP; FLJ33877; RC-GAP72; RCGAP72; p73; p73RhoGAP ; FILGAP

Immunogen:

Human ARHGAP24.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A