

ARHGAP26 Antibody



PACO07824

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGAP26

Uniprot

Q9UNA1

Synonyms:

Rho GTPase activating protein 26; ARHGAP26; FLJ42530; GRAF; KIAA0621; OPHN1L; OPHN1L1 ;

Immunogen:

Human ARHGAP26.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A