

ARHGAP5 Antibody



PACO07826

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGAP5

Uniprot

Q13017

Synonyms:

Rho GTPase activating protein 5; ARHGAP5; GFI2; RhoGAP5; p190-B ;

Immunogen:

Human ARHGAP5.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C,
Avoid freeze / thaw cycles.

Product Images

N/A

N/A