

# ARHGDIB Antibody



PACO07828

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARHGDIB

**Uniprot**

P52566

**Synonyms:**

Rho GDP dissociation inhibitor (GDI) beta; ARHGDIB; D4; GDIA2; GDID4; LYGDI; Ly-GDI; RAP1GN1; RhoGDI2 ;

**Immunogen:**

Human ARHGDIB.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A