

ARHGEF16 Antibody



PACO07833

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ARHGEF16

Uniprot

Q5VV41

Synonyms:

Rho guanine exchange factor (GEF) 16; ARHGEF16; GEF16; NBR ;

Immunogen:

Human ARHGEF16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A