

ARL15 Antibody



PACO07840

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ARL15

Uniprot

Q9NXU5

Synonyms:

ADP-ribosylation factor-like 15; ARL15; ARFRP2; FLJ20051 ;

Immunogen:

Human ARL15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A