

ARL3 Antibody



PACO07843

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ARL3

Uniprot

P36405

Synonyms:

ADP-ribosylation factor-like 3; ARL3; ARFL3 ;

Immunogen:

Human ARL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A