

ARL5B Antibody



PACO07846

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ARL5B

Uniprot

Q96KC2

Synonyms:

ADP-ribosylation factor-like 5B; ARL5B; ARL8 ;

Immunogen:

Human ARL5B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A