

ARL6IP1 Antibody



PACO07848

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARL6IP1

Uniprot

Q15041

Synonyms:

ADP-ribosylation factor-like 6 interacting protein 1; ARL6IP1; AIP1; ARL6IP; ARMER; KIAA0069 ;

Immunogen:

Human ARL6IP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A