

ARL6IP6 Antibody



PACO07851

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARL6IP6

Uniprot

Q8N6S5

Synonyms:

ADP-ribosylation-like factor 6 interacting protein 6; ARL6IP6; MGC33864; PFAAP1 ;

Immunogen:

Human ARL6IP6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A