

ARMC8 Antibody



PACO07856

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ARMC8

Uniprot

Q8IUR7

Synonyms:

armadillo repeat containing 8; ARMC8; HSPC056; MGC10058; MGC4880; S863-2 ;

Immunogen:

Human ARMC8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A