

# ARPC3 Antibody



PACO07861

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARPC3

**Uniprot**

O15145

**Synonyms:**

actin related protein 2/3 complex, subunit 3, 21kDa; ARPC3; ARC21; p21-Arc ;

**Immunogen:**

Human ARPC3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A