

ARR3 Antibody



PACO07864

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ARR3

Uniprot

P36575

Synonyms:

arrestin 3, retinal (X-arrestin); ARR3; ARRX ;

Immunogen:

Human ARR3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A