

ASB8 Antibody



PACO07875

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ASB8

Uniprot

Q9H765

Synonyms:

ankyrin repeat and SOCS box-containing 8; ASB8; FLJ21255; MGC5540; PP14212 ;

Immunogen:

Human ASB8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A