

ASGR1 Antibody



PACO07882

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ASGR1

Uniprot

P07306

Synonyms:

asialoglycoprotein receptor 1; ASGR1; ASGPR; CLEC4H1; Hs.12056 ;

Immunogen:

Human ASGR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A