

PACO07885

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ASMTL

Uniprot

O95671

Synonyms:

acetylserotonin O-methyltransferase-like; ASMTL; ASMTLX; ASMTLY; ASTML ;

Immunogen:

Human ASMTL

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A