

ASPRV1 Antibody



PACO07890

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ASPRV1

Uniprot

Q53RT3

Synonyms:

aspartic peptidase, retroviral-like 1; ASPRV1; MUNO; SASP; SASPase; Taps ;

Immunogen:

Human ASPRV1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A