

ATF6B Antibody



PACO07899

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ATF6B

Uniprot

Q99941

Synonyms:

activating transcription factor 6 beta; ATF6B; CREB-RP; CREBL1; FLJ10066; G13 ;

Immunogen:

Human ATF6B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A