

ATL1 Antibody



PACO07909

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATL1

Uniprot

Q8WXF7

Synonyms:

atlastin GTPase 1; ATL1; AD-FSP; FSP1; GBP3; SPG3; SPG3A; atlastin1 ;

Immunogen:

Human ATL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A