

ATP2C1 Antibody



PACO07921

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATP2C1

Uniprot

P98194

Synonyms:

ATPase, Ca++ transporting, type 2C, member 1; ATP2C1; ATP2C1A; BCPM; HHD; KIAA1347; PMR1; SPCA1; hSPCA1 ;

Immunogen:

Human ATP2C1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A