

ATP5H Antibody



PACO07926

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ATP5H

Uniprot

O75947

Synonyms:

ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit d; ATP5H; ATP5JD; ATPQ ;

Immunogen:

Human ATP5H.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A