

ATP5I Antibody



PACO07927

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP5I

Uniprot

P56385

Synonyms:

ATP synthase, H+ transporting, mitochondrial F0 complex, subunit E; ATP5I; ATP5K; MGC12532 ;

Immunogen:

Human ATP5I.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A