

# ATP6AP1 Antibody



PACO07931

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ATP6AP1

**Uniprot**

Q15904

**Synonyms:**

ATPase, H+ transporting, lysosomal accessory protein 1; ATP6AP1;16A; ATP6IP1; ATP6S1; Ac45; CF2; MGC129781; VATPS1; XAP-3; XAP3 ;

**Immunogen:**

Human ATP6AP1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A