

ATP6V1C1 Antibody



PACO07936

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP6V1C1

Uniprot

P21283

Synonyms:

ATPase, H+ transporting, lysosomal 42kDa, V1 subunit C1; ATP6V1C1; ATP6C; ATP6D; FLJ20057; VATC; Vma5 ;

Immunogen:

Human ATP6V1C1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A