

ATP6V1D Antibody



PACO07938

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP6V1D

Uniprot

Q9Y5K8

Synonyms:

ATPase, H+ transporting, lysosomal 34kDa, V1 subunit D; ATP6V1D; ATP6M; VATD; VMA8 ;

Immunogen:

Human ATP6V1D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A