

# ATP6V1F Antibody



PACO07941

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ATP6V1F

**Uniprot**

Q16864

**Synonyms:**

ATPase, H+ transporting, lysosomal 14kDa, V1 subunit F; ATP6V1F; ATP6S14; MGC117321; MGC126037; MGC126038; VATF; Vma7 ;

**Immunogen:**

Human ATP6V1F.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A