

ATP6V1G1 Antibody



PACO07942

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP6V1G1

Uniprot

O75348

Synonyms:

ATPase, H+ transporting, lysosomal 13kDa, V1 subunit G1; ATP6V1G1; ATP6G; ATP6G1; ATP6GL; ATP6J; DKFZp547P234; Vma10 ;

Immunogen:

Human ATP6V1G1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A