

ATP6V1G3 Antibody



PACO07943

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ATP6V1G3

Uniprot

Q96LB4

Synonyms:

ATPase, H+ transporting, lysosomal 13kDa, V1 subunit G3; ATP6V1G3; ATP6G3; MGC119810; MGC119813; Vma10 ;

Immunogen:

Human ATP6V1G3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A