

ATP8A1 Antibody



PACO07944

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATP8A1

Uniprot

Q9Y2Q0

Synonyms:

ATPase, aminophospholipid transporter (APLT), class I, type 8A, member 1; ATP8A1; ATPASEII; ATP1A; ATP2; MGC130042; MGC130043; MGC26327 ;

Immunogen:

Human ATP8A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A