

ATP9A Antibody



PACO07945

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATP9A

Uniprot

O75110

Synonyms:

ATPase, class II, type 9A; ATP9A; ATP1IA; KIAA0611 ;

Immunogen:

Human ATP9A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A