

# ATPAF1 Antibody



PACO07946

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ATPAF1

**Uniprot**

Q5TC12

**Synonyms:**

ATP synthase mitochondrial F1 complex assembly factor 1; ATPAF1; ATP11; ATP11p; FLJ22351; MGC88060 ;

**Immunogen:**

Human ATPAF1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A