

ATPIF1 Antibody



PACO07948

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATPIF1

Uniprot

Q9UII2

Synonyms:

ATPase inhibitory factor 1; ATPIF1; ATPi; ATPiP; IP; MGC1167; MGC8898 ;

Immunogen:

Human ATPIF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A