

ATR Antibody



PACO07949

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ATR

Uniprot

Q13535

Synonyms:

ataxia telangiectasia and Rad3 related; ATR; FRP1; MEC1; SCKL; SCKL1 ;

Immunogen:

Human ATR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A