

# B3GALNT2 Antibody



PACO07960

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

B3GALNT2

**Uniprot**

Q8NCR0

**Synonyms:**

beta-1,3-N-acetylgalactosaminyltransferase 2; B3GALNT2; B3GalNAc-T2; MGC39558 ;

**Immunogen:**

Human B3GALNT2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A