

B3GALT6 Antibody



PACO07962

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

B3GALT6

Uniprot

Q96L58

Synonyms:

UDP-Gal: betaGal beta 1,3-galactosyltransferase polypeptide 6; B3GALT6; beta3GalT6 ;

Immunogen:

Human B3GALT6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A