

B4GALT3 Antibody



PACO07967

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

B4GALT3

Uniprot

O60512

Synonyms:

UDP-Gal: betaGlcNAc beta 1,4- galactosyltransferase, polypeptide 3; B4GALT3; beta4Gal-T3 ;

Immunogen:

Human B4GALT3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A