

BCAP31 Antibody



PACO07984

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

BCAP31

Uniprot

P51572

Synonyms:

B-cell receptor-associated protein 31; BCAP31;6C6-AG; BAP31; CDM; DXS1357E ;

Immunogen:

Human BCAP31.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A